



POLAR OH1

Getting Started Guide

COMPATIBLE WITH



Model: 2L


support.polar.com/en/OH1

MANUFACTURED BY


POLAR ELECTRO OY
PROFESSORINTIE 5
FI-90440 KEMPELE
TEL +358 8 5202 100

www.polar.com


1.




2.




3.




4.



5.



6.



This guide will help you get started with Polar OH1. For more help and product information see support.polar.com/en/OH1. Polar OH1 can be used as a training device that records and stores data or as a heart rate sensor with a Polar wrist unit.

Charge

Place the sensor into the USB adapter with the lens facing up so that the contacts on the sensor and USB adapter meet (picture 1). Plug the USB adapter into a computer USB port or USB power wall outlet (picture 2). Do not charge the device when the charging contacts are wet.

Use with Polar Beat mobile app

Download Polar Beat onto your mobile device. Turn OH1 on by pressing the button. Pair the sensor from Beat settings and follow instructions.

Use with Polar wrist unit

Download Polar Flow app onto your mobile device. Turn OH1 on by pressing the button. Flow app discovers the sensor automatically and guides you through the setup. If you want to use a computer, install Polar FlowSync from flow.polar.com/start, run it, plug the sensor into a USB port and follow instructions. Remember to pair the sensor with your wrist unit. See support.polar.com for device-specific pairing instructions.

Wear on your arm

Place the sensor into the armband holder with the lens facing up (picture 3). Wear the armband so that the sensor is on the underside of the armband firmly against your skin (picture 4). Wear the OH1 around your lower or upper arm, not around your wrist (picture 5).

Wear with swimming goggles

Place the sensor in one of the included swimming holders, clip the holder onto your swimming goggles strap (picture 6). The lens should touch the bare skin of your temple.

Start training

You can record your training session with your wrist unit, Polar Beat app or with just the OH1 sensor alone. If you want to use the sensor alone, turn OH1 on by pressing the button until the lights are switched on and press the button twice to start the session. When you're done training press and hold the button until the lights are switched off.



By default, OH1 sends your HR signal also via ANT+. You can turn it off from the Polar Beat settings, if you want.

Maintenance

See the armband care label for washing instructions. Treat the sensor with care.

Materials

Sensor: ABS, ABS+GF, PMMA, SUS 316 (Stainless steel).



This product is not a toy. Do not allow children or pets to play with this product. This product contains small components that can be a choking hazard.

The radio equipment operates 2.402 - 2.480 GHz ISM frequency band(s) and 2.0 mW maximum power.