AGU[®] ADVANCED
GROWING
UP

COMPRESSOR
NEBULIZER

BALLOON

EN INSTRUCTION
MANUAL

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TABLE OF CONTENTS

1	Introduction	3
2	Safety precautions	4
3	Product identification. Durable period	6
4	Technology of adjusting the valve	7
5	Assembling your nebulizer kit	8
6	Operating your nebulizer	9
7	Using the AC Adapter	10
8	Cleaning procedures	11
9	Changing the air filter	12
10	Troubleshooting	13
11	Specifications	14
12	Technical data	15

1 INTRODUCTION

Thank you for purchasing an **AGU N5** Compressor Nebulizer. With proper care and use, your nebulizer will provide you with many years of reliable treatments.

This unit operates on standard AC power. Treatments are delivered quickly, safely and conveniently making this unit ideal for all ages. We encourage you to thoroughly read this guidebook to learn about the features of your nebulizer.

Your Compressor nebulizer **AGU N5** relates to medical equipment and should be used under the supervision of a licensed physician and / or a respiratory therapist.

Together with your physician and / or respiratory therapist, you can feel comfortable and confident knowing that you are obtaining the most effective inhalation treatments for your respiratory condition.

NOTE: Your nebulizer is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy. Please consult with your physician and / or pharmacist to determine if your prescription medication is approved for use with this nebulizer. For type, dose, and regime of medication follow the instructions of your doctor or respiratory therapist.

This device fulfils the provision of the EC directive 93 / 42 / EEC (Medical Device Directive) and the European Standard EN 13544-1:2007+A1:2009 Respiratory therapy equipment – Part 1: Nebulizing systems and their components.

Please read this manual carefully before use and be sure to keep this manual.

CAUTIONS

Please use general safety precautions when operating your nebulizer. This unit should be used only for its intended purpose as described in this guidebook and with medications only under the supervision and instruction of your physician. Do not use the device in anesthetic or ventilator breathing circuits.

2 SAFETY PRECAUTIONS

PRODUCT CAUTIONS

Read the following before using

- To avoid electrical shock: keep unit away from water.
- Do not immerse the unit in liquid.
- Do not use while bathing.
- Do not reach for a unit that has fallen into water immediately unplug the unit.
- Do not use the unit if it has any damaged parts (including plug), if it has been submerged in water or dropped. Promptly send the unit for examination and repair.
- The unit should not be used where flammable gas, oxygen or aerosol spray products are being used.
- Keep the air vents open. Do not place the unit on a soft surface where the openings can be blocked.
- If the medication cup is empty, do not attempt to operate the unit.
- If any abnormality occurs, discontinue to use until the unit has been examined and repaired.
- The unit should not be left unattended while plugged in.
- Do not tilt or shake the unit when in operation.
- Disconnect the unit from the electrical outlet before cleaning, filling and after each use.
- Do not use attachments unless recommended by the manufacturer.
- Do not plug or unplug the AC adapter into the electrical outlet with wet hands.
- Unplug the AC adapter from the electrical outlet after using the device.
- Do not disassemble or attempt to repair the Compressor Unit.
- Do not use the device in anesthetic or ventilator breathing circuits.

OPERATING CAUTIONS

- Close adult supervision is highly recommended when the unit is used by children or disabled person.
- Keep your eyes away from the output of medication mist.

- The maximum capacity of the medication cup is 5 ml and should not be overfilled.
- Do not use this unit while operating a vehicle.
- If any discomfort or abnormality occurs, stop using the unit immediately.
- Do not use the device if the air tube is bent.

STORAGE CAUTIONS

- Do not store the unit in direct sunlight, high temperature or humidity.
- Keep the unit out of reach of small children.
- Always keep the unit unplugged while not in use

CLEANING CAUTIONS

- Check air filter, mouthpiece, nebulizer and any other optional components before each use. Dirty or worn parts should be replaced.
- Do not immerse the unit in water. It may damage the unit.
- Disconnect the unit from the electrical outlet before cleaning.
- Clean all necessary parts after each use as instructed in this guidebook.
- Always dispose of any remaining medication in the medication cup after each use. Use fresh medication each time you use the device.
- Do not store the air tube with moisture or medication remaining in the air tube. This could result in infection as a result of bacteria.


MEDICAL DISCLAIMER

This manual and product are not meant to be a substitute for advice provided by your doctor or other medical professionals. Don't use the information contained herein or this product for diagnosing or treating a health problem or prescribing any medication. If you have or suspect that you have a medical problem, promptly consult your doctor.

3

PRODUCT IDENTIFICATION. DURABLE PERIOD

Nebulizer kit


Nozzle

Air filter


Child mask


Angled mouthpiece

Air tube


AC adapter

Air filter slot

Air output

AC adapter jack


Power switch

Durable periods are as follows, provided the product is used to nebulize 2 ml of medication 2 times a day for 10 minutes each time at room temperature (25 °C).

Durable period may vary depending on usage environment.

- Main unit 5 years
- Nebulizer kit 1 year
- Air tube, mouthpiece 1 year
- Air filter 60 days
- Child mask 5 years.

4 TECHNOLOGY OF ADJUSTING THE VALVE

The proprietary adjustable valve is able to deliver medications of different viscosity level according to every user's conditions and needs at ease. Our VA technology allows users to adjust different levels of nebulization rate ranging from 0.08 – 0.2 ml/min at consistent particle size.


The nebulization rates can be adjusted by the user in a very easy way without exchanging parts. Higher nebulization rate/fully open is for higher viscosity medications and higher breathing capacity user while lower nebulization rate with closed valve will be more appropriate for kids/infants with lower breathing capacity.

5

ASSEMBLING YOUR NEBULIZER KIT

Follow the cleaning instructions in this guidebook under «Cleaning technique» prior to using your nebulizer for the first time or after it has been stored for an extended period of time.

REMEMBER: Always unplug the compressor and make sure the power switch is turned to the «OFF» position before cleaning, assembling and before or after each use.


1. Place the compressor on a flat, stable surface within reach.
2. Gently twist and pull straight up on the lid of the nebulizer to separate into two parts (medication cup and cover).
3. Make sure that the nozzle is properly installed inside the medication cup. The stem inside the medication cup inserts into the tube of the nozzle.
4. Reassemble the nebulizer by carefully twisting the medication cup and cover together. Make sure that the two parts fit securely.

6

OPERATING YOUR NEBULIZER

The nebulizer is operable at up to a 45° angle. If the angle is greater than 45°, no aerosol will be generated.

1. Add the prescribed amount of medication to the medication cup.
2. Attach one end of the air tube connector to the air output.
3. Carefully attach the opposite end of the air tubing connector to the stem at the base of the nebulizer kit.
4. Attach the angled mouthpiece or mask to the top of the nebulizer.

NOTE: A 30-minute interval is recommended after each use.

Operating your nebulizer

After every use:




- Unplug the unit from the power source.
- Allow the unit to completely cool.
- Carefully detach the air tubing from the nebulizer and pour out any remaining medication.
- Follow the cleaning procedures provided in this guidebook.

7 USING THE AC ADAPTER

AC Adapter




($\varnothing 3.8/\varnothing 1.35$)

1. Check that the power switch is in the «OFF» position.
2. Connect the AC adapter with the AC adapter jack in the back of the unit (with required voltage and current indicated near the AC adapter jack).
3. Plug the AC adapter into the socket.
4. Use only the authorized AC adapter with this nebulizer.
5. Press the power switch to the «ON» position.

8 CLEANING PROCEDURES

Rinsing technique (performed after each treatment or before first use).

1. Make sure that the power-switch has been turned to the «OFF» position and the unit has been disconnected from the power source.
2. Disconnect the air tube from the nebulizer device.
3. Gently twist and pull up the cover of the nebulizer kit to open and separate.
4. Rinse the nebulizer kit and components with hot tap water.
5. Dry with clean towels or completely air dry.
6. Reassemble the nebulizer kit.

NOTE: For the first time cleaning or after the unit has been stored for an extended period of time, thoroughly clean all components, including the air tube. The nebulizer kit is dishwasher safe.


Cleaning the compressor

Wipe the compressor daily using a soft cloth.

NOTE: Any other form of cleaning or cleaning agents may damage the finish of the unit.

9 CHANGING THE AIR FILTER

It is important to change the air filter approximately when the air filter turns gray.

1. Remove the air filter cover by gently pulling forward.
2. Discard the gray filter.
3. Replace with a new, clean air filter.
4. Securely reattach the air filter cover to the unit.

NOTE: Air filters cannot be cleaned or washed. Only **AGU N5** air filters can be used. Do not substitute alternate material such as cotton. Do not operate without an air filter.

10 TROUBLESHOOTING

If any abnormality occurs during use, please check and correct the following.

1. Unit does not operate when power switch is pressed:
 - Check the AC connection to the outlet.
2. No misting or low rate of misting:
 - Check that there is medication in the nebulizer cup.
 - Check the main unit if there is any physical damage.
 - Check the position of the nozzle inside the nebulizer.
 - Make sure that air tube and other components are properly attached.
 - Check the air filter and replace if necessary.

Protection against electric shock:

- Class II equipment.

Type BF applied parts:

- Mouthpiece and mask.

Protection against harmful ingress of water and particulate matter:

- IP21

Degree of safety in the presence of flammable anesthetics or oxygen:

- No AP / APG (not suitable for use in the presence of flammable anesthetics or oxygen).

11 SPECIFICATIONS

Type	Compressor nebulizer AGU N5
Model	AGU N5  NB60
AC adapter	DC12 V, 1 A
Power	AC 100~240 V, 50~60 Hz
Power consumption	≤ 12 W
Sound level	≤ 55 dBA (1 meter away from AGU N5)
Compressor pressure range	> 16 psi (110 kPa)
Operating pressure range	> 8 psi (55 kPa)
Operating flow range	> 2.0 lpm
Operating temperature range	10 °C to 40 °C
Operating humidity range	≤ 90% relative maximum humidity
Operating atmospheric pressure range	700–1060 hPa
Storage temperature range	-20 °C to 60 °C
Storage humidity range	≤ 90 % RH
Dimension (L x W x H)	134.7 x 88.3 x 54 mm
Weight	286.5 g (without accessories)
Medication capacity	5 ml (cc)
Particle size (MMAD)	< 3.0 μm
Average nebulization rate	Fully open valve > 0.2 ml / min (0.9% saline solution) Closed Valve > 0.08 ml / min (0.9% saline solution)
Standard accessories	Nebulizer kit, air tube, mouthpiece, filters (5 pcs), child mask

Subject to technical modification without prior notice.


Performance may vary with drugs such as suspensions or high viscosity, see drug supplier's data sheet for further details.

12 TECHNICAL DATA

According to the Standard EN 13544-1:2009+A1:2009 «Respiratory therapy equipments Part 1: Nebulizing systems and their components», Annex CC.3 using the multistage cascade impactor method to measure the particle size.

Particle Size	MMAD $\leq 3 \mu\text{m}$ (MMAD = Mass Median Aerodynamic Diameter)
Medication Cup Capacity	5 ml max.

Particle Size Distribution

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC EMISSIONS

The **AGU N5** is intended for use in the electromagnetic environment specified below.

The customer or the user of the **AGU N5** should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment—guidance
RF emissions CISPR 11	Group 1	The AGU N5 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The AGU N5 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Compliance	

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC IMMUNITY

The **AGU N5** is intended for use in the electromagnetic environment specified below.

The customer or the user of the **AGU N5** should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+ 6 kV contact + 8 kV air	+ 6 kV contact + 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient / burst IEC 61000-4-4	+ 2 kV for power supply lines + 1 kV for input / output lines	+ 2 kV for power supply lines Not applicable	Main power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	+ 1 kV line(s) to line(s) + 2 kV line(s) to earth	+ 1 kV differential mode Not applicable	Main power quality should be that of a typical commercial or hospital environment
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec.	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AGU N5 requires continued operation during power mains interruptions, it is recommended that the AGU N5 be powered from an uninterruptible power supply or a battery
Power frequency (50 / 60 Hz) magnetic field IEC 61000-4-8	3 A / m	3 A / m	The AGU N5 power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

NOTE: UT is the a.c. mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC IMMUNITY

The **AGU N5** is intended for use in the electromagnetic environment specified below. The customer or the user of the **AGU N5** should assure that is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the AGU N5 including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V / m 80 MHz to 2,5 GHz	3 V / m	

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations or radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **AGU N5** is used exceed the applicable RF compliance level above, the **AGU N5** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **AGU N5**.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

RECOMMENDED SEPARATION DISTANCE between portable and mobile RF communications equipment and the AGU N5.

The **AGU N5** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **AGU N5** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **AGU N5** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARNING: The symbol on this product means that it's an electronic product and following the European directive 2012/19/EU the electronic products have to be disposed on your local recycling centre for safe treatment.

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Manufacturer: Rossmax Swiss GmbH (factory Rossmax (Shanghai) Incorporation Ltd. in China), SWITZERLAND

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Date of manufacturing is on the device in the YYMMDD format:

YY – year of manufacturing; MM – month of manufacturing;

DD – day of manufacturing. Made in China.

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Make sure that children could not use the appliance without supervision, because some small parts can get into the respiratory tract or the digestive system.

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The set completeness and operation of the device are tested in the presence of the buyer.

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ATTENTION!

This warranty is valid only if documents are filled in properly.

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Applications to service and repair center.

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AGU[®] is the registered trademark by Montex Swiss AG,