Avicenum **LYMPHO 360**

(EN) Instruction manual

Avicenum LYMPHO 360 Arm Sleeve

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Avicenum LYMPHO 360 Armstrumpf

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(GB)

INSTRUCTION MANUAL

Dear customer, the medical compression arm sleeve that you hold in your hands is a useful aid in the treatment of circulatory and lymphatic disorders. You can maximize the benefits of this product through proper use and maintenance.

The Avicenum LYMPHO 360 Arm Sleeve is a class II compression sleeve (II.KT, CCL2). It was made by circular knitting technology. The materials used in our compression arm sleeves are kind to the skin and do not contain latex. The compression arm sleeve works by exerting a precisely defined pressure on the arm, within the range 23–32 mmHg. The proper fit is necessary to achieve the desired effect. Use the compression sleeve when engaged in physical activity; compression products should be removed for longer rest periods (unless otherwise indicated by your doctor). **The Avicenum LYMPHO 360 Arm Sleeve** is a medical aid intended for use both at home and in medical care facilities.

Product name:

Avicenum LYMPHO 360 Arm Sleeve

Product variants:

- Gloveless, hem, class II, size S normal – XL long
- Glove, hem, class II, size S normal – XL long
- Glove, shoulder support, class II, size S normal – XL long

User profile:

The Avicenum LYMPHO 360 Arm Sleeve is a medical aid intended for use both at home and in medical care facilities.

Intended use of the medical aid:

This medical aid is intended for compression therapy of circulatory disorders and various types of oedema in the arm.

Indications:

- lymphoedema
- oedemas caused by parasitic disease
- malignant and benign tumours of lymphatic vessels (angiomas and sarcomas)
- lipoedema, lipo-lymphoedema
- massive localised lymphoedema accompanying extreme obesity
- inflammation of the lymph nodes (microbial, post-traumatic, postoperative)
- malignant lymphomas
- chronic venous insufficiency
- general oedemas among women

If you are using a compression arm sleeve for the first time, consult your doctor to determine the proper compression class and have regular check-ups during treatment.

Contraindications:

- cardiac oedema
- peripheral ischemic artery disease
- arterial hypertension
- heart rhythm disorders
- ischemic heart disease
- progressive systemic sclerosis (scleroderma)
- chronic polyarthritis
- sudeck's atrophy (Complex regional pain syndrome)
- hypersensitivity to materials used

In view of known contraindications your doctor should evaluate possible risks and weigh them against the therapeutic benefits of use when considering application of the compression arm sleeve.

Target group:

Patients with lymphoedema, primary or secondary oedema, and patients with the above indications. The compression arm sleeve is intended for use on healthy, intact skin. This compression arm sleeve is suitable for all age groups and both sexes. There are no restrictions for use during pregnancy.

Side effects:

Side effects may occur, especially if the arm sleeve does not fit properly or if it is worn when at rest. There may be bruising, impaired blood flow due to pressure on blood vessels, or impaired sensitivity due to excessive nerve compression. Sensitive individuals may experience skin irritation (redness, itching, or blistering) in places where the compression arm sleeve comes into contact with the skin.


Fitting your compression arm sleeve:

The proper fit is necessary to achieve the desired effect. Use the size chart to determine the correct size by measuring limb circumference at the points indicated.

Note:

Any arm oedema should be examined by your doctor. If this is the first time you experience oedema, consult your doctor before using the product. Do not use the product if you have not determined the cause of the oedema. Always consult your doctor before using this product in combination with other medical aids (such as bandages or orthoses).

- use a tape measure
- measure first thing in the morning—limbs may swell over the course of the day
- ask someone to help you
- measurements may change over the course of treatment – always check your size before using a new compression arm sleeve

SIZE	Circumference of the hand in the measured points			
	a	c	d	g
S	up to 18	16 - 19	22 - 26	24 - 30
M	up to 20	19 - 22	26 - 30	30 - 36
L	up to 23	22 - 25	30 - 34	36 - 42
XL	up to 25	25 - 28	34 - 38	42 - 48
Length	normal		long	
C - G	up to 45		over 45	

- a** palm circumference
c wrist circumference
d forearm circumference
g underarm circumference
C – G length

If the measured values do not correspond to any of the sizes in the table (e.g. point c falls under size S, but point g is size L), visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression arm sleeve.

If your measurements fall between sizes, choose the larger size.

If you have difficulties putting on the compression arm sleeve, visit a specialised medical supplies store and request a demonstration of the correct method of putting it on. There are also aids available to help you put it on.

If the hem does not stay put on the arm, clean the silicone studs with an alcohol-based solution. Shave the area where the silicone comes into contact with the skin.

If the compression arm sleeve slips down, tighten the shoulder support.

Correct fit:

- You should feel noticeable but not unpleasant pressure over the entire length of the arm sleeve
- The compression arm sleeve should not slip or roll over at the elbow
- The hem of the compression arm sleeve should not be too tight
- The shoulder support of the compression arm sleeve should not be too tight

With a correctly selected size:

- you should feel noticeable but not unpleasant pressure over the entire length of the product
- the product does not slide nor roll down below the knee or on the instep when walking
- the hem does not choke

Information about compression class, size, and manufacturer can be found on the tag sewn at the upper edge of the compression arm sleeve. Should the tag be removed in the course of use, sizing information is also knit in at the bottom of the sleeve.

Putting on the compression arm sleeve:

Pull the compression sleeve onto the arm gradually; do not tug on the upper edge alone. It is best to put on the compression arm sleeve in the morning before any oedema occurs. When pulling on the compression arm sleeve, keep it away from sharp objects (e.g. jewellery) or long nails. The compression sleeve must be stretched evenly, without folding, to ensure effective pressure distribution. If you use any body cosmetics, allow 15 minutes after application before putting on the sleeve. In the event of visible mechanical damage to the product, stop wearing it and replace it with a new one.

Care and maintenance:

This compression arm sleeve is a medical aid intended for repeated use by a single person. The product will maintain its properties and safety for 6 months from the first use under these conditions:

- the compression arm sleeve may be machine washed (max. temperature 30 °C, gentle cycle) with a soap solution or a special detergent for compression medical aids, in a protective cover (washing bag)
- do not use fabric softener
- do not spin-dry, press the wet compression sleeve dry between two towels
- lay flat to dry away from a direct heat source (radiator, direct sunlight)
- do not iron the compression sleeve
- maximum number of washing cycles is 120
- store in a dark, dry place in the original packaging
- do not alter the compression arm sleeve in any way, do not repair areas damaged by use, do not trim the edges or hem
- to refresh the adhesive properties of the silicone studs along the hem, wipe clean with an alcohol solution
- pay special attention to the hook-and-loop fastener (for compression arm sleeves with shoulder support):
 - contact of the hook-and-loop fastener with the knit area may cause damage
 - contact of the hook-and-loop fastener with other materials (fabric, textiles, dirt and grime) may damage it and cause it to stop working in whole or in part
- keep the compression arm sleeve away from organic solvents and bleaching agents

This medical aid may be used safely for 60 months from the date of manufacture. This date is part of the product batch number expressed in the format YY/MM-1234567, where YY is the year and MM the month of manufacture.

Care symbols:

Do not use fabric softeners

When not to use:

Stop using the product immediately if any of the following conditions apply:

- There are signs of mechanical damage to the sleeve – replace with a new compression arm sleeve
- You notice a deterioration in your health that may be associated with the use of the compression arm sleeve – consult your doctor in such case

Material:**Avicenum LYMPHO 360 Arm Sleeve**

- Gloveless, hem: 60 % PAD nylon, 40 % elastane LYCRA®, silicone
- Glove, hem: 60 % PAD nylon, 40 % elastane LYCRA®, silicone
- Glove, shoulder support: 60 % PAD nylon, 40 % elastane LYCRA®

This medicinal aid does not contain an active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, or tissues or cells of animal origin or derivatives thereof.

Disposal:

The compression arm sleeve can be disposed of with normal household waste or recycled as a textile. Although the product is intended for use on healthy, intact skin, it may become contaminated with body fluids during use. In that case, dispose of the compression arm sleeve as infectious material. The material used is non-toxic unless ignited. Never incinerate used compression medical aids.

Any serious adverse event occurring in connection with the product in question should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located:

Czech Republic	Státní ústav pro kontrolu léčiv, Šrobárova 48, 100 41 Prague 10
Slovak Republic	Štátny ústav pre kontrolu liečiv, Kvetná 1024/11, 821 08 Bratislava
Poland	Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Das Bundesinstitut für Arzneimittel und Medizinprodukte, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kiev, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

List of abbreviations:

KT, CCL – compression class

mmHg – millimetres of mercury column (unit of blood pressure measurement on the leg)

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We hereby certify that

Aries A.S.

is certified as a License Partner for the region of the Czech Republic, and that it is certified to a region quality control program.

The company License Partner is authorized to use the Aries logo and has the right to distribute its respective articles with the Certified Quality Seal.

OEKO-TEX®
CERTIFICATE

The company

Aries, A.S.
Studenec 309
674 01 Studenec u Prahy, CZECH REPUBLIC

A production facility according to EN ISO 9001:2015 is in use (including the year 2015) since September 1st 2015.

for the following articles

Woolen, cotton and synthetic goods (textiles, accessories, certified good products, and wool fabric) (including the year 2015) in the Czech Republic, and good products and wool fabric (including the year 2015) in the Czech Republic, and good products and wool fabric (including the year 2015) in the Czech Republic.

Products of the respective type according to EN ISO 9001:2015, according to a regional quality control program for the respective region (including the year 2015) in the Czech Republic.

As certified partner, the License Partner is authorized to use the Aries logo and has the right to distribute its respective articles with the Certified Quality Seal.

Customer of the certificate, who has received a certificate according to EN ISO 9001:2015 is obliged to use the Aries logo and to cooperate with products that are not in the scope of the certificate.

The certificate is valid until **31.03.2018**.

NO. 14.02.001

Signature: _____
Date: _____

Aries
MEDICAL PRODUCTS

ARIES, a.s., Studenec 309
512 33, MADE IN CZECH REPUBLIC
ISO 9001:2015

www.avicenum.eu

SD 11-005 (R04)