AT53064

**ORTEZA KRĘGOSŁUPA INSTRUKCJA UŻYTKOWANIA
SPINAL ORTHOSIS USER MANUAL
BENUTZERHANDBUCH FÜR WIRBELSÄULENORTHESSEN
NÁVOD K POUŽITÍ ORTÉZA PÁTEŘE
NÁVOD NA POUŽITIE ORTÉZA CHRBTICE
MANUEL D'UTILISATION DE L'ORTHÈSE RACHIDIENNE
SPINALE ORTHESE GEBRUIKERSHANDLEIDING
MANUAL DE USUARIO DE LA ÓRTESIS ESPINAL
SPINALE ORTHESE GEBRUIKERSHANDLEIDING
ANVÄNDARHANDBOK FÖR SPINALORTOS**

MD

PL

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ENG

GENERAL WARNING

Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product until you have read and understood these instructions. If you do not understand the warnings, cautions and recommendations, please contact your healthcare professional or dealer to avoid damage to the product.

NOTE: Inspect all parts for damage caused in transit. If such damage is noticed, DO NOT use the product. Contact the manufacturer for more information.

APPLICATION

The spinal orthosis is used to stabilise and strengthen the spine in the sacral, lumbar and thoracic sections.

INDICATIONS

- lumbosacral pain syndromes
- osteoporosis of the spine
- increased muscle tension
- lumbar discopathy
- spinal instability

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling. Injuries, inflammation, open wounds, abrasions, degeneration that prevent comfortable use of the product.

FEATURES

The brace consists of two parts: a high lumbar and sacral part and a Velcro-fastened dorsal part, which allows the appropriate size to be adjusted to the patient's height.

SIZE (To choose the right size, measure your waist circumference)

M: 70 - 85 cm

L: 80 - 95 cm

XL: 90 - 105 cm

XXL: 100 - 115 cm

XXXL: 110 - 125 cm

XXXXL :120-135 cm

Height with dorsal section: 55 cm

Height of lumbar and sacral section: 26 cm

SCOPE OF SUPPLY

Spinal brace (2 parts), instruction manual

NOTES:

In the event of a product-related 'serious incident' that directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) the death of a patient, user or other person; or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

this 'serious incident' should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions or other distressing, unclear symptoms related to the use of the medical device, you should stop using the device and consult a healthcare professional.

ATTENTION: It is forbidden to use the product other than for its intended purpose

ATTENTION: The manufacturer is not responsible for damage caused by negligent cleaning and maintenance, or resulting from failure to follow the instructions in this manual.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this orthosis is intended (see section on intended use of the stabilisator in these instructions). The product can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/sizes/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

USE

1. At the very beginning, the high part of the orthosis should be attached to the lower part with Velcro.
2. Make sure the orthosis is in the correct position, fitting the orthosis according to the illustrations below.
3. The brace is put on by fastening the main strap with Velcro, the strap should be pressed tightly enough.
4. Then grasp both additional drawstrings at the same time and fasten on the main belt, adjusting the correct tightness.
5. The additional straps located on the high part of the orthosis should be placed on the shoulders, crossed at the back and fastened to the main strap at the front, the covers should be placed under the armpits.
6. Finally, it is important to make sure once again that the device fits properly on the spine, is worn at the right height and that the straps are tight enough to provide stabilisation.

NOTE: The device can be used using only the lower part (omitting the step of attaching the upper part) as long as it is in accordance with the doctor's recommendations.

CLEANING AND MAINTENANCE

The product is recommended to be washed by hand with a mild detergent, at a temperature of up to 40°C. Do not wash in a washing machine. Do not use strong detergents, bleach or rough objects for cleaning. Dry flat laid out in the shade, do not expose to direct sunlight. Do not tumble dry. Do not dry clean. Do not iron.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

HOW TO DISPOSE OF THE PRODUCT

Once the orthosis has been decommissioned, the medical device can be disposed of as normal municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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