ENGLISH

Electric toothbrush

Description

- Travel case
- Battery compartment lid
- Body
- Attachments
- Base
- Timer
- Cup

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Before using the unit, read these instructions carefully. Keep these instructions for future reference.

Use the unit for intended purposes only, as it is stated in this manual.

Mishandling the unit can lead to its breakage and cause harm to the user or damage to his/her property.

SAFETY MEASURES

- Do not wash the toothbrush attachments in a dishwashing machine.
- If necessary, consult your dentist before using this unit.
- You should consult your dentist before using this unit if your child has undergone oral cavity or gum surgery in the last 2 months.
- If after using the toothbrush there is a profuse bleeding and if it does not stop during a week, consult your dentist.
- If your child has a heart pacemaker or any other device implanted, consult your doctor or the manufacturer of the implanted device before using the unit.
- This unit is intended for tooth and gum brushing. Use the unit only for the intended purposes.
- If you feel any discomfort or pain during tooth brushing, stop using the unit and consult your doctor.
- A toothbrush is a personal hygiene item that should only be used by one individual.
- Use only the attachments supplied or the attachments recommended by the manufacturer.
- Wash your toothbrush thoroughly after every usage, especially if your toothpaste contains peroxide, cooking soda and other bicarbonates (which are often used in whitening toothpastes). It prevents possible plastic cracking.
- Do not immerse the body of the toothbrush into water or other liquids. It may lead to unit damage and breakage.

THE UNIT IS INTENDED FOR HOUSEHOLD USAGE ONLY

Base (5) assembling

- Insert the ledges on the upper part of the base in the openings on the lower part.
- Use the base for brush (3), attachments (4), timer (6) and cup (7) storage.
- You can mount the base (5) on the wall.

USAGE

- Before the first usage wash the attachment (4) with warm water during 30 seconds (pic. 1).
- Set the attachment (4) on the body (3) until the clamp clicks.
- Set the battery compartment lid (2) clamp to the position ● (pic. 2) and take the battery compartment lid (2) off.
- Insert 2xAAA batteries (supplied), set the lid (2) back to its place, set the clamp to the position ▲.
- Switch the unit on by pressing the button «I».
- Put the rotating head of the toothbrush to your tooth, hold it in such position for a few seconds.
- Moving the head of the toothbrush from one tooth to another clean the outer surface of your upper teeth (30 seconds), then inner surface of your upper teeth (30 seconds). Clean outer and inner surfaces of your lower teeth (for 30 seconds each).

Note: Recommended total time for tooth brushing is 2 minutes. Use the timer (6) for measuring time.

- Clean the chewing surface of your teeth.
- Switch the unit off by pressing the button «0».

CLEANING AND CARE

- Wash the attachment (4) under warm water jet after every usage.
- Make extra cleaning of the attachment not less than once a week. For this purpose take the attachment off and wash it thoroughly in and outside.
- Do not immerse the body of the toothbrush (3) into water or other liquids. It may lead to unit damage and breakage.
- IMPORTANT: Replace the attachment (4) in case of deformation and wear of bristles (at least than every 3 months). Only high quality bristles ensure effective and thorough cleaning of your teeth.**
- Use only the attachments for toothbrushes (interchangeable attachments WX-2251 ST).

STORAGE

- Clean the unit before taking it away for storage.
- Use the travel case (1) for long storage and transportation of the unit.
- Keep the unit away from children in a dry cool place.

DELIVERY SET

Toothbrush with attachment – 1pc.
Interchangeable attachment – 1pc.
Timer – 1 pc.
Cup – 1 pc.
Base – 1 pc.
Travel case – 1 pc.
«AAA» batteries - 2 pcs.
Instruction manual – 1 pc.

TECHNICAL SPECIFICATIONs

Power supply: 2 «AAA» batteries, 1,5 V

The manufacturer preserves the right to change the specifications of the unit without a preliminary notification.

Unit operating life is 1 year

Guarantee

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC)

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© ООО ГОЛДЕР-ЭЛЕКТРОНИКС, 2012
© GOLDBER-ELECTRONICS LLC, 2012