SB-316

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ENGLISH

Please read this instruction manual carefully before using the appliance.

Save the instruction manual as you may require it in future.

This device is designed to determine the body weight by weight, which affects it.

The device is intended for use only for domestic purposes in accordance with this Instruction Manual.

After unit transportation or storage at low temperature it is necessary to keep it for at least 4 hours at room temperature before switching on.

TECHNICAL SPECIFICATIONS

Maximum weight: 150 kg

Minimum weight: 2,5 kg

Graduation: 100 g

Power: 2 x 1,5V AAA

Service lifetime - 3 years

Warranty period - 12 months

IMPORTANT SAFEGUARDS

Attention! Incorrect installation of batteries may result in an explosion. Batteries should only be replaced with the same or equivalent type.

- Read this manual carefully before operating the scales to avoid damage during use. Incorrect handling can lead to damage to the product, cause material damage or cause damage to the health of the user.
- The scales are not intended for use by persons (including children) with reduced physical, sensory or mental abilities, or lack of experience or knowledge, unless they are under the control or instructed about the use of the device by the person responsible for their safety.
- Scales purchased in the cold season, in order to avoid failure, must be kept for at least four hours at room temperature.
- Do not use the device for purposes not specified in the instruction manual.
- Avoid getting water on the scales to prevent the unit from breaking down.
- Place the scales on a flat surface. Do not stand on the edge of the scales. Stand only on the central part of the platform.
- It is forbidden to stand on the scale to more than one person at a time.
- Do not stand on the scale if the surface is wet. This can lead to slipping.
- Do not use scales near sources of electromagnetic radiation.
- Do not overload or subject the balance to impactloads.
- Do not jump on the scales.
- Do not allow the packaging material to fall into the hands of children (risk of suffocation).
- Before switching on, make sure the batteries are in good condition and properly installed.

ENGLISH

- Use the batteries specified in the instructions. Do not use batteries with visible signs of oxidation.
- Do not charge batteries, as this may cause an explosion.
- Do not dispose of batteries with other litter. Pass them to special recycling centers.
- Use the scales only for household purposes in accordance with this Instruction Manual. Scales are not intended for industrial or commercial use. For malfunctions caused by the use of the product for purposes beyond the scope of personal domestic needs (ie for industrial or commercial purposes) is not covered by the guarantee.
- The manufacturer reserves the right to make minor changes to the design without additional notice products that do not fundamentally affect its safety, performance and functionality.

INSTALLING AND REPLACING OF BATTERIES

Attention!

Danger of explosion! Do not expose batteries to high temperatures, direct sunlight or fire. Do not dispose batteries in a fire. The battery contains chemicals, so it must be disposed of in accordance with the applicable disposal regulations. To prevent leakage of the electrolyte, which could result in bodily injury, property damage or damage to the balance, observe the following rules:

- Install 2 AAA batteries correctly, observing the polarity indicated at the battery installation location (+ and - signs).
- Remove the batteries if the scales aren't used for a long time.

OPERATION

Attention! To obtain accurate data, install the scales only on a firm, flat, horizontal surface and take readings in one place. Do not place the scales on the carpet. Before the first use, remove the electrical insulation tape from the lithium battery.

If "Lo" is displayed on the LED display, open the battery compartment and replace it with a new one.

If the permissible weight (150 kg) is exceeded, "EEEE" will be displayed on the display. Immediately step off the scales when this message appears.

For weighing, select the necessary weight unit (the weight unit is located on the rear panel) and stand on the scales. After a few seconds, your weight will be showed on the display. The scales switch off automatically after you leave them.

ADDITIONAL FUNCTIONS

In order to measure the percentage of water, fat, muscle and bone mass in the body, you need to do the following:

1. Press the "SET" button. The settings menu will open.
2. Enter the user number (from 0 to 9) using the "▲" and "▼" buttons, and then confirm the entries by pressing the "SET" button.
3. Use the "▲" and "▼" buttons to select your gender, then confirm the entries by pressing the "SET" button.
4. Use the "▲" and "▼" buttons to enter your height, then confirm the entries by pressing the "SET" button.
5. Use the "▲" and "▼" buttons to enter your age, then confirm the entries by pressing the "SET" button.

Display indication	Meaning
kg/lb/st	Weight (kg, lb, st)
P0-9	User number
	Male
	Female
cm/ft/in	Height (cm, ft, in)
Age	Age
FAT	Fat
TBW	Water
MUS	Muscles
BONE	Bone
%	Percentage

6. After 2 seconds, the display will show "0.0" and you will be able to take readings.

7. Stand on the scales. Make sure that your feet are dry, and the heels are evenly positioned on the metal plates-electrodes. The display alternately displays your weight, percentage of fat, water, muscle and bone mass, as well as the number of kilocalories that you need to consume daily for the normal functioning of your body.

8. Remember your number. The next time you take the water, fat, muscle and bone mass readings, you just need to choose the howling number, wait until the display shows "0.0", and you can take the reading.

DISPLAY INDICATION

Below is a table with which you can find out what form you are in.

Age	Male		Female		Fat quantity	Display indication
	Fat	Water	Fat	Water		
≤30	4.0-20.5	66.0-54.7	4.0-15.5	66.0-58.1	Lack	-
	20.6-25.0	54.6-57.6	15.6-20.0	58.0-55.0	Norm	0
	25.1-30.5	51.5-47.8	20.1-24.5	54.9-51.9	Excess	+
	30.6-50.0	47.7-34.4	24.6-50.0	51.8-34.4	Overweight	++
>30	4.0-25.0	66.0-51.6	4.0-19.5	66.0-55.4	Lack	-
	25.1-30.0	51.5-48.1	19.6-24.0	55.3-52.3	Norm	0
	30.1-35.0	48.0-44.7	24.1-28.5	52.2-49.2	Excess	+
	30.1-50.0	44.6-34.4	28.6-50.0	49.1-34.4	Overweight	++

ENGLISH

CAUTION: If the display shows an error message, wait until the device turns off. Please enable it again to restart the operation. If the error message does not disappear, remove the batteries from the battery compartment for 10 seconds, then insert them back and turn on the device.

Display indication	Meaning
8888	Select user and personal settings.
----	Incorrect position of the scales. Move the scales to a straight, hard surface.
Lo	The batteries are exhausted. It is necessary to replace them with new ones.
Err	The percentage of fat is outside the minimum or maximum value of the measuring scales.
EEEE	The scales are overloaded. Weight exceeds the maximum permissible weight for weighing. Get off the platform to avoid damaging the scales.

CLEANING, MAINTENANCE AND STORAGE

TRANSPORTATION. Transportation of the device is necessary by any kind of covered transport with the application of cargo securing rules ensuring preservation of the presentation of the product and / or packaging and its further safe operation. **DO NOT** expose the device to shock loads during handling operations.

CLEANING. It is recommended to wipe the scales with a damp cloth. Do not use abrasive cleaners, detergents containing powder, alcohol, gasoline, or various chemicals. Do not stand on the scale in shoes to avoid damage to the scales.

STORAGE. The device must be stored indoors, in conditions that presume the preservation of the presentation of the product and its further safe operation.

WARRANTY TERMS

When purchasing the product, ask for it to be checked in your presence. Without the having a warranty card or when it is not properly filled in, the quality claims are not accepted, and no warranty repair is performed. The warranty card should be presented at any address to the service center during the whole warranty period. The product is accepted for service in a complete set only. The warranty period is calculated from the moment of sale to the buyer. Keep the documents confirming the date of purchase of the product (cash receipt). The condition of free warranty service of your product is its correct operation, which does not go beyond personal domestic needs, in accordance with the requirements of the operating instructions for the product, the absence of mechanical damages and the consequences of careless handling of the product.

The product is accepted for warranty service in its pure form. The warranty covers all manufacturing and construction defects (except those listed in the "Warranty does not apply" section), identified during the warranty period. During this period defective parts, except consumables, are subject to free replacement in the warranty service center.

THE GUARANTEE DOES NOT COVER:

1. Defects caused by force majeure.
2. Damage to the product caused by the use of the product for purposes that go beyond personal domestic needs (ie industrial or commercial purposes).
3. Consumables and accessories (sections, etc.).
4. Defects caused by overloading, improper operation, penetration of liquids, dust of insects, ingress of foreign objects into the product.
5. Products that have been repaired outside authorized service centers or warranty workshops.
6. Damage resulting from changes in the design of the product by the user or unqualified replacement of his components.
7. Mechanical damages, including due to negligent handling, improper transportation and storage, falling of the product.
8. Violation of the operating instructions.
9. Incorrect installation of mains voltage (if required).
10. Making technical changes.
11. Damage through the fault of animals (including rodents and insects).

Production month: refer to the information on the gift box.
For warranty service, please contact your nearest service center.

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