INFRARED FOOT MASSAGER**MODEL: DH 77L****IMPORTANT**

- read all instruction manual before the first operating and save it
- never immerse the appliance, cord and plug into water
- never use the appliance with wet hands
- never use the appliance on a damp floor
- never use the appliance with an extension lead
- always unplug the appliance before assembling/disassembling
- always unplug the appliance before cleaning
- before the first operating make sure that the voltage in your wall outlet is the same as on the appliance's rating label marked
- never use the damaged appliance or appliance with damaged cord/plug
- never use the appliance with the damaged wall outlet
- any examination, repair or adjustment must be done in an authorized Service Center
- never use the appliance near or on the hot surfaces
- never let the cord hang over edge of the table or touch any hot surfaces
- use only attachment recommended by the manufacturer
- always unplug the appliance when not in use
- never let the children play with the appliance
- give a high attention when the appliance is being used by or near the children
- use only the appliance for it's direct purpose
- never expose the appliance to fall down or to be kicked
- never use caustic or abrasive cleaners

2

*IM_hydromassager_DH77L.indd*²
Process Black

- This appliance has a heating function. People insensitive to heat must be specially careful during use.

- If you have any concerns regarding your health, consult your doctor before using the massager.

- Individuals with pacemakers should consult a doctor before using the massager.

- In case of pregnancy, diabetes or illness, consult your doctor before using the massager.

- We recommend to use the massager for not longer than 20 minutes at a time. Extensive use could cause the product to overheat. Should this occur, discontinue use and allow the massager to cool before operating.

- If you feel pain in muscles or become pain for a prolonged period of time, stop using the massager and consult your doctor. Persistent pain could be a symptom of a more serious condition.

- Never use any massager or infrared heat on open wounds, discolored areas, or any area of the body that is swollen, burned, inflamed or where skin eruptions or sores are present.

- Massage and infrared heat should be pleasant and comfortable. Should pain or discomfort result, stop using and consult your doctor.

SPECIAL WARNINGS

- Never immerse the massager in water or other fluids.

- Always check the temperature of the water before placing your feet into the massager.

- Do not stand on the massager. Sit down and place your feet in the massager.

- Always unplug the massager when it is not in use and when filling, emptying, cleaning or moving the massager.

- Ensure that your hands are dry when operating and switching.

- Do not leave the massager unattended while it is operating.

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IM_hydromassager_DH77L.indd

4
Process Black**ENGLISH**

Heating function will only keep the temperature of water placed inside the massager! The massager will not heat cold water! For a warm water massage, please fill the massager with warm water!

OPERATING

- ensure that the massager is unplugged before filling with water - do not fill water over the MAX-mark
- select the function by pressing the corresponding button on the massager or on the remote control - the corresponding LED will illuminate.
- 1 - Bubbles/Heat
- 2 - Massage/Heat
- 3 - Infrared
- 4 - OFF (all functions are switched off)
- sit down and place your feet in the spa
- enjoy a deep penetrating massage by sliding feet back and forward over the reflex rollers

PEDICURE CENTRE

- 1 - select desired attachment
 - pumice stone**
to smooth and soften dry skin on heels and soles of feet
 - brush**
to clean and exfoliate skin and for soft massage
 - massage balls**
enjoy a relaxing massage from your heels to the balls of your feet
- 2 - Place the attachment on the centre post
- 3 - push down with your foot to activate the attachment
 - to empty water - unplug the massager and tip the massager allowing water to spill out along the drain spout

ATTENTION!

Never pour water over the control panel and the remote control

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22.04.2004, 12:10

ENGLISH**CLEANING**

- Unplug the massager and allow it to cool before cleaning - clean only with a soft, damp cloth.
 - Never allow water or any other liquid to contact the control panel and the remote control.
 - After cleaning, wipe surfaces with a dry cloth.
 - Rollers can be removed for cleaning - place thumb and forefinger around roller and pull out.
 - Never use abrasives, detergents or cleaning agents.
 - Periodically straightens the cord if it becomes twisted.
- STORAGE**
- Unplug the massager from the outlet and allow it to cool.
 - Store in its box in a clean and dry place.
 - To avoid breakage, do not wrap the power cord around the massager.
 - Do not hang the massager by the power cord.
- This appliance is designed for household use only!**

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