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ATTENTION! SAVE THIS OPERATING MANUAL FOR FURTHER USE!

Dear Buyer! Thank you for purchasing Planta products. The myostimulator EMS-750 is a versatile innovative portable body massager using an electrical stimulation method, ideally designed for convenient use.

The advantage of this kit for myostimulator of PLANTA EMS-750 mice is the presence of additional accessories (belt, elbow/knuckle, slippers) in the kit, 10 pieces of replaceable reusable (from 70 to 100 uses) silicon electrodes (ellipse shape), 2 pieces of which are spare.

This myostimulator PLANTA EMS-750 compact, you can always take it with you.

You can always use the device both for one user and for two users at the same time, also using the accessories in the kit.

The electrical stimulation method is based on the procedure of electrical stimulation of muscles and nerves, carried out by transmitting pulse current with specified characteristics from the myostimulator to the human body through electrodes. Electrical stimulation as a separate direction in medicine and cosmetology has been formed for many centuries. Back in ancient times, people used the electric action of amber and discharges of electric fish to restore the muscle activity of damaged tissues and eliminate diseases associated with nerves.

Today, electrical stimulation methods are widely used in rehabilitation medicine, sports, as well as cosmetology and weight loss. The principle of operation of the electrical stimulator is based on the simulation of autogenic electrical impulses, which are conducted through the skin to nerve or muscle fibers using electrodes, acting on sensitive nerve endings.

Advantages PLANTA EMS-750:

- Combining 2 stimulation techniques TENS and EMS
- 6 massage zones, 4 kinds of massage, 2 massage programs
- 2 channels – independently adjustable, with the ability to use the maximum number of electrodes simultaneously in various configurations with additional accessories (belt, armrest/knuckle, slippers)
- ability to use on two users at the same time
- 3 kinds of using: 8 electrodes or 4 electrodes and accessory or 2 accessories at the same time
- Large and accurate LCD-display
- Intuitive and simple using of device
- Possibility of setting massage time
- Extensive set of accessories – SLIPPERS, WAIST, KNEE/ELBOW, 5 pairs of electrodes
- Power saving function
- Compact size

TENS & EMS

Thanks to the combination of two effective methods at once – myostimulation (EMS) and electroneurostimulation (TENS), the massager relaxes perfectly, relieves fatigue and eliminates discomfort.

Electroneurostimulation (TENS) is the skin-acting electrostimulation of nerves and nerve endings. Weak electrical impulses (no more than 200 ms), affecting sensitive nerve endings, block pain impulses, eliminating the sensation of stiffness, and also increase the production of “endorphins”

– joy hormones that improve mood, have an analgesic and calming effect.

The neurostimulation method for 2017 is considered the most effective way to suppress chronic pain syndrome. Pain is just a signal that is transmitted through nerves from a problem area to the brain. Pain is information to the brain that there is a malfunction in the body and it needs to be eliminated. But this “just as information” can make a person’s life unbearable, especially when it comes to chronic pain.

Electromyostimulation (EMS) – muscle electrostimulation, which gives two effects at once: relaxation and stimulation. In the first case, muscle relaxation occurs, tension and fatigue are relieved, which is very important with increased physical activity. In the second case, the muscles are toned, their elasticity is restored, endurance is increased. Under the influence of impulse current, muscles begin to contract, as a result of which oxygen-enriched blood is more actively poured into them, lymph outflow occurs and muscle tone is improved.

Thanks to two separately adjustable channels, the device allows the pulse intensity to be adjusted independently of each other, for example, to cover the body on both sides or evenly stimulate large areas of tissues. Individual adjustment of the intensity of each channel makes it possible to simultaneously use the device in different parts of the body, due to which time savings are achieved compared to sequential use.

Using effect

Regular use of this device contributes to:

- removing stress and fatigue, relieving emotional stress;
- maintaining the tone of the mice;
- getting rid of stiffness and severity in muscles with increased physical activity and after training, or vice versa with a sedentary lifestyle and sedentary work;
- elimination of congestion and improvement of blood circulation;
- improvement of quality of sleep and general well-being, replenishment of life forces.

! SAFETY REQUIREMENTS

READ ALL INSTRUCTIONS BEFORE USING!

When using this device, always observe the basic safety rules:

- Use the instrument solely for the purposes for which it was designed and only in the manner described in this application manual.
- This device is powered by batteries and mains. Therefore, it is subject to the same safety rules as any other electrical appliances. Connect the power supply unit of the device only to the mains with the appropriate voltage.
- Ensure that children do not use the device as a toy. The use of the device by children and persons with limited physical, sensory or mental abilities, or persons who lack the knowledge and experience necessary for the use of the device, is allowed only if the person responsible for the safety of the above mentioned category of users instructs them about the safe use of the device and controls their actions, and this category of users is aware of the risks associated with the use of the device.
- Ensure that all instrument parts, electrodes, wiring, accessories, power supply are in good condition before use. Use only the original accessories included. If a defect or malfunction is

detected, immediately turn off the device and stop using it.

- Only authorized service specialists can repair the device.
- Never open or disassemble the device yourself.
- Individual components of this device cannot be disassembled, repaired or reinstalled.
- Do not use the device while sleeping, driving machines and equipment or other work, or under the influence of alcohol or other substances.
- The device is intended for external use only.
- Do not apply electrodes to the anterior surface of the chest (the area between the ribs and the sternum), primarily in the area of both large thoracic muscles. This may increase the risk of atrial flicker.
- Do not apply electrodes to the bones of the skull, in the area of the eyes, mouth, pharynx or larynx, in the area of the neck/carotid arteries, in the area of the genitals.
- Do not remove electrodes when the device is on and the wires are inserted, always disconnect the device or the corresponding channel before moving or removing the electrodes.
- Do not confuse the electrode cable with contacts with headphones or other devices, do not connect the electrodes to other devices. Do not modify the electrodes (for example, trimming them), this will lead to an increase in current density.
- Pre-treatment of the skin with oily cream or ointment is not recommended, this will greatly increase the wear of the electrodes and can lead to current spikes.
- Do not use pointed objects (handles, toothpicks, etc.) to press the buttons of the device.
- Protect the sample from falls.
- Keep the device in a clean, dry, cool place, protected from direct exposure to sunlight.
- In no case immerse the device and its accessories in water or other liquids!
- To avoid device breakage, do not use the device and do not touch the connecting wires of the electrodes and the power cord with wet hands.
- Do not use a myostimulator near short and microwave devices (for example, mobile phones), as this can lead to current jumps.
- Never use this device in conjunction with electromedical equipment, electronic life support systems such as mechanical ventilation devices, or portable medical devices such as an electrocardiograph. This device can cause malfunctions in the operation of the above-mentioned medical equipment.
- Ensure that metal objects, such as belt buckles or chains, do not contact the electrodes during stimulation.
- If decorations or metal objects are present in the application area, they must be removed before the device can be used in order to avoid burns.

Counterindication:

If you have any serious chronic diseases, you take any drugs, or you have any contraindications, then before using the device you should consult your doctor. Immediately stop using the device if you feel discomfort or weakness.

It is forbidden to use the device to those people who have:


- skin sensitivity disorders (e.g., reduced pain threshold), any skin diseases and inflammation, allergic reactions on the skin, as well as open wounds, bruises and postoperative scarring during healing

- cardiovascular disorders, heart attack, stroke
- tumors, oncology, inflammation of internal organs and tissues
- infectious diseases, sepsis, elevated body temperature
- hypertensive and other vascular diseases
- gastrointestinal diseases, gallstones or kidneys
- blood clotting disorders, any diseases of veins and lymph nodes, thrombophlebitis and tendency to thromboembolic diseases
- pacemaker, insulin pump, implants, prostheses, pins or helix
- for diseases occurring with seizures (e.g. epilepsy)
- Diabetes

It is forbidden to use the device for menstruation, pregnancy and after surgery!

⚠ ATTENTION! THIS INSTRUMENT MODEL NOT A MEDICAL PRODUCT!

DETAILS & COMPONENTS

1. “ON/OFF” button.

2/18 Adjustment of “CHANNEL 1” and “CHANNEL 2” intensity.

Rotation clockwise increases the intensity of the action, rotation counterclockwise – reduces the intensity. Total of 20 intensity levels are available.

Buttons No. 3-5, 10-12 – selection of massage area:

“ARM” “JOINT,” “FEET,” “SHOULDER,” “BACK,” “LEGS”

Pressing the button once enables massage mode for the corresponding area. The selected mode is displayed on the display.

Buttons № 4-9 – selection of massage type:

“INTENSIVE,” “ACUPUNCTURE,” “VIBRO MASSAGE,” “THUMP”.

Pressing the button once includes the corresponding type of massage.

The selected mode is displayed on the display.

13. “WEIGHT CONTROL” massage program selection button

14. “FITNESS” massage program selection button

15. Socket for connection of wires with electrodes and accessories “CHANNEL 1”

16. Socket for connection of wires with “CHANNEL 2” electrodes and accessories

17. “+” and “-” buttons to adjust massage time

19. Power cord hole

20. Battery compartment (no batteries included)

ADDITIONAL ACCESSORIES

1.

6.

3.

4.

5.

8.

2.

7.

1. 1 x PLANTA myostimulator EMS-750

2. 10 x silicon gel electrodes (reusable, ellipse shape)


3. 2 x wires per 4 electrodes

4. 2 x wires for accessories

- 5. 1 x BELT
- 6. 1 x PAIR OF SLIPPERS
- 7. 1 x BOW BAR/BOW BAR
- 8. 1 x power adapter 6V

NOTE: The manufacturer reserves the right to change parts of the product without prior notice for different schedule lines without affecting the basic technical parameters of the product. This can entail a change in the weight and dimensions of the product, but not more than $\pm 5-10\%$.

DISPLAY

1. Intensity indicator "CHANNEL 1."
Displays the selected intensity level for Channel 1 – 1 to 20
2. Intensity indicator "CHANNEL 2"
Displays the selected Channel 2 intensity level – 1 to 20
3. "INTENSIVE" massage indicator
4. "ACUPUNCTURE" massage indicator
5. "VIBRO MASSAGE" massage indicator
6. Indicator for selection of massage type "THUMP"
7. Timer
Shows the time remaining until the massage session ends
8. Indication of massage mode for the selected area of electrodes application – "HANDS" "JOINT," "FEET," "SHOULDER," "BACK," "LEG". When "HANDS" mode is selected, indication of this area and signature "ARM"; "SHOULDER" - "shoulder"; "BACK" - "waist"; "LEG" - "leg";
9. Indicator of instrument operation (flashes when used)


10. Indication of "WEIGHT CONTROL" program selection

11. Indication of selection of "FITNESS" program

Note! In "WEIGHT CHECK" and "FITNESS" programs the modes change, the currently active mode will be highlighted on the display accordingly.

FOR 2 PEOPLE AT THE SAME TIME

The device EMS-750 is a model with 2 independently controlled channels. Thanks to this, you have the opportunity to connect accessories to the electric stimulator in different configurations for one or two people at the same time.


FIRST USING

- Remove the device and all accessories from the box. Remove all packaging materials and labels.
- Carefully check the device and all accessories for defects and damages. If you find defects and damage, do not use the device and contact the seller.
- During the first use, it is recommended to test the device sitting on a chair or lying down to get used to the current.
- The operating intensity should be adjusted according to individual sensations.


MANUAL

1. Connect the electrodes or accessory in various configurations (see "Combinations") to the connection wires as shown in the figures below. Tightly insert the connector wire plug into the electrode connector so that the bare contact metal is not visible. Otherwise, the contacts of the device will not be connected, the device will not work. To connect accessories, wires with buttons are used, which snap on the connectors located on the accessories.


Electrodes and accessories have various connecting wires (see Section- Equipment and Accessories)

Pic. 1 Connection electrodes

Pic.2 Connection waist & knee/elbow

Pic.3. Connecting the slippers

2. Connect wires to the device & accessories as shown in the figure below. Before connecting the connection wires to the device, make sure that the device is turned off. Hold the insulated parts of the connecting wire, insert the connecting wire plug into the connector of the corresponding channel on the lower part of the device housing. Verify that the wires are connected correctly. The myostimulator has two channels. You can use one channel or two channels together by connecting connection wires to one or two connectors. Using two channels gives the user the advantage of stimulating different parts of the body at the same time, and also makes it possible to use the device at once by 2 users at the same time.


USING OF STICKERS- ELECTRODES

Remove protective from electrodes surface.


Attention!

Do not discard protective films because they will need to be glued back to the electrodes at the end of the procedure. Detailed electrode care recommendations are listed in the "Care and Maintenance" section!

Apply the electrode with the adhesive side to a clean dry and smooth area of the skin and press it tightly to the skin. Ensure that the electrodes are aligned and securely fixed to the skin. If you select an electrode location on the body, follow the electrode location rules in the Electrode Location section. Do not turn on the device until the electrodes are properly attached to the skin.

3. To enable the device, insert the network adapter into the connector on the device housing and connect the myostimulator to the network. Then press "ON/OFF button" on the main part of the device housing, as shown in the figure below, at that the digital LCD display of the device will light up.

Only swap the devices - first extinguished screen, then lit.

4. Select the application and one of the massage types, or one of the WEIGHT CONTROL or FITNESS programs. You can set different massage modes separately for channel 1 and channel 2.

5. Set the required massage time from 5 to 25 minutes by pressing the "+" or "-" buttons sequentially. The "+" button increases the massage time by 1 minute, the "-" button reduces the set massage time by 1 minute.

The massage time is displayed at the top of the LCD and defaults to 20 minutes.

The myostimulator will automatically shut down after the set time. You can also complete the massage at any time by clicking on the "ON/OFF".

6. Set the desired intensity level from 1 to 20 for the appropriate channel.

Intensity levels of "CHANNEL 1" and "CHANNEL 2" are displayed on the display on both sides in the form of scales.

7. To pause, press the ON/OFF button once, and the electrodes will pause. Press the button once to continue the massage "ON/OFF," the electrodes will start working again.

8. After the massage is complete, turn off the myostimulator by pressing the "ON/OFF" button.


Disconnect the connecting wires from the myostimulator and electrodes and remove the electrodes from the skin. To remove the electrodes, lift the electrode by the corner and carefully disconnect it from the skin. Glue the protective film to the adhesive side of the electrode.

COMBINATIONS OF ELECTRODE AND ACCESSORY APPLICATIONS

The user of the device may use any possible combination of electrodes and accessories. Electrodes and accessories can be applied to one person or at once to two users.

Below you will find the main applications of electrodes and accessories, but you can also use other combinations on your choice.

EXAMPLES OF COMBINATION ADDITIONAL PARTS

USING OF ADDITIONAL ACCESSOIRES

Sequentially switch on and connect accessories, as well as switch off device control, refer to section "Operating Manual"

To use accessories (WAIST, SLIPPERS, KNEE/ELBOW), use certain wires (see section "Packing and additional accessories")

Each accessory has fasteners, connecting which you can connect the accessory.

Slippers

To use, it is necessary to dress slippers, connect them to the device, select one of the modes of operation "INTENSIVE," "ACUPUNCTURE," "VIBRO MASSAGE" or "THUMP" the area of application "FEET" and set a comfortable intensity.

Application Effect:

Slippers are a special accessory that effectively relieves painful tension in the muscles, relieves the decay of the legs, improves coordination of the movement of the toes, reduces the convulsions of the foot and shin mice.

It is recommended for treatment and prevention to athletes, people whose profession is associated with a long stay on their feet, as well as before and after intensive physical training.

Belt

For use, it is necessary to dress the belt, connect to the device, select one of the modes of operation "INTENSIVE," "ACUPUNCTURE," "VIBRO MASSAGE" or "THUMP" the area of application "BACK" and set a comfortable intensity.

Application Effect:

Belt is a special accessory with a wide spector of use (as a myostimulator for the press, as well as on the lumbar spine, hip joint, hip). It improves coordination of lumbar movement, increases circulation in the lumbar region, strengthens the lumbar-sacral spine and pumps the muscles of the press.

KNEE/ELBOW

For use it is necessary to put on an elbow/knee piece kneecap (knee/elbow), to choose one of operating modes "INTENSIVE", " VIBRO MASSAGE" or "THUMP", «, a scope of "HAND", "JOINT" or "LEG" and to establish to VIBRO comfortable intensity.

Application Effect:

Knee/elbow - accessory significantly expands the possibilities of exposure to electric therapy devices, allows you to most effectively exert a therapeutic effect on the problem zone. By contracting smooth muscles in the walls of the vessels activates cellular respiration, metabolism. It is recommended for athletes, people whose profession is associated with a long stay on their feet, as well as before and after intensive physical training. It improves the coordination of the movement of the legs and hands, increases the circulation of the ulnar and knee joints, reduces "crunches" and "clicks" during movement.

SELECTION OF LABEL ELECTRODES LOCATION

Proper placement of electrodes is important for the successful application of electrical stimulation. When selecting the electrode location, observe the following instructions:


Distance between electrodes:

The greater the distance between the electrodes, the greater the area of the tissue to be stimulated. This applies to flat and concave areas of the body. However, as the distance between the electrodes increases, the tissue stimulation force decreases. This means that by choosing a larger distance between the electrodes, a larger area is stimulated, but with less force. To increase stimulation, it is necessary to increase the intensity of the pulse. When selecting the distance between the electrodes, follow the instructions:

- the distance between the electrodes should be not less than 7 cm, if the distance is less than 7 cm, this will lead to too strong stimulation of surface structures.
- recommended distance between electrodes from 7 to 15 cm
- the distance between electrodes greater than 15 cm will lead to very weak stimulation of the surface with a large area.

Arrangement of electrodes depending on direction of muscle fibers

The efficiency of electrical stimulation in all modes largely depends on the correct choice of the location of the electrodes on the body. Ideally, the electrodes should be placed on the motor points of the muscle or nerve. The motor point of the nerve is the area where the nerve is most superficially located under the skin and accessible by current. The motor point of the muscle is the place corresponding to the place where the motor nerve enters the muscle. The choice of current flow direction must be consistent with the direction of muscle fibers.


- If superficial muscles are to be reached, the electrodes should be arranged parallel to the direction of the fibers (A - B/C - D),
- If deep layers of fabrics are to be reached, the electrodes must be arranged across the fiber direction. The latter can be achieved by a cruciform (transverse) arrangement of electrodes, for example, A - D/B - C. With electrostimulation of nerves and nerve endings (TENS) the electrodes of one channel should be installed so that the point at which there is discomfort is located between the electrodes, or one electrode would be located directly at the point at which there is discomfort, and the second electrode - at a distance of at least 7 cm from the first. The electrodes of the remaining channels can be used to simultaneously work out other points at which there is discomfort or, together with the electrodes of the first channel, used to block the entire area where there is discomfort (on the opposite side). Here again, a cruciform arrangement makes sense.

ELECTRODE LAYOUT (EMS/TENS)

To use the EMS function, select the modes "WEIGHT CONTROL," "FITNESS," "ACUPUNCTURE," "HAND" "JOINT" "FEET" "SHOULDER" "BACK," "LEG". To use the TENS function, select the programs "INTENSIVE", "ACUPUNCTURE", "VIBRO MASSAGE".

Be obligatory to place the electrodes in the correct cruciform!

Example:


ELECTRODE OVERLAY OPTIONS (EMS/TENS)

The following table shows the arrangement of electrodes or accessories.

	Programs	tens/ems	Electrode's location (pic)	Resistance 0.5K	Frequency	Pulse Width
				(VP)	(Hz)	(US)
1.	ARM	ems	Pic.12,13,14,15,16,17 - Electrodes; KNEE/ELBOW PAD	102V	1Hz	280uS
				102V	2.48Hz	280uS
				102V	4.96Hz	280uS
				102V	15.5uS	280uS
				78V	49uS	280uS
				102V	19uS	280uS
2.	JOINT	ems	Pic.1,12,14,16,17,23,24 25,28 Electrodes; KNEE/ELBOW PAD or SLIPPERS	102V	1Hz	280uS
				102V	2.48Hz	280uS
				102V	4.96Hz	280uS
				102V	19uS	180uS
				78V	49uS	180uS
				102V	15.5uS	50-200uS
3.	FEET	ems	Pic.28 - Electrodes; KNEE/ELBOW PAD & SLIPPERS	102V	4.9Hz	280uS
				102V	1Hz	280uS
				102V	2.48Hz	280uS
				102V	19uS	180uS
				102V	15.5uS	50-200uS
4.	SHOULDER	ems	Pic.12,13,14,15,16,17 - Electrodes; KNEE/ELBOW PAD	102V	1Hz	280uS
				102V	2.48Hz	280uS
				102V	4.96Hz	180uS
				102V	19Hz	180uS
				78V	49Hz	180uS
				102V	15.5Hz	50-200uS
5.	BACK	ems	Pic.3,4,5,6,7,8,9,10,11 Electrodes; WAIST	102V	1Hz	280uS
				102V	2.48Hz	280uS
				102V	4.96Hz	280uS
				102V	19uS	180uS
				78V	49uS	180uS
				102V	15.5uS	50-200uS
6.	LEG	ems	Pic.22,23,24,25,26,27,28 Electrodes; KNEE/ELBOW PAD or SLIPPERS	102V	1Hz	280uS
				102V	2.48Hz	280uS
				102V	4.96Hz	180uS


7.	INTENSIVE	tens	This program can be used on any part of the body	78V	49Hz	180uS
8.	ACUPUNCTURE	tens	This program can be used on any part of the body	102V	15.5Hz	50-200uS
9.	VIBRO MASSAGE	tens	This program can be used on any part of the body	102V	19Hz	180uS
10.	THUMP	ems	This program can be used on any part of the body	102V	1Hz	280uS
11.	WEIGHT CONTROL	ems	Рис.18,19,20,21 KNEE/ELBOW PAD or WAIST	102V	4.96Hz	180uS
				102V	15.5Hz	50-200uS
				102V	19Hz	180uS
				78V	49Hz	180uS
12.	FITNESS	ems	Рис.18,19,20,21 KNEE/ELBOW PAD or WAIST	102V	1Hz	280uS
				102V	4.96Hz	180uS
				102V	2.48Hz	180uS

Muscles of neck zone

This method of using electrodes works out deep muscles of the cervical-collar zone. Facilitates the study and preparation of mice for physical activity. Used to relieve fatigue, maintain the tone of mice. Select the massage area and type of massage "INTENSIVE," set the timer for 15 minutes, or "ACUPUNCTURE" for 5-10 minutes. Perform massage 2-3 times a week.


Muscles of back

This method of using electrodes improves blood circulation, processes and tones the muscles of the back and lower back.

Select the massage area and type of massage "INTENSIVE," set the timer for 15 minutes, or "ACUPUNCTURE" for 10 minutes. Perform massage 2-3 times a week.

Muscles of arm & shoulder

This method of using electrodes improves blood circulation, processes and tones the muscles of the hands and forearms. Facilitates the study and preparation of mice for physical activity. It is also used to relieve fatigue, maintain the tone of mice.

Select the massage area and type of massage "INTENSIVE," set the timer for 15 minutes or "ACUPUNCTURE" or "VIBRO MASSAGE" for 5-10 minutes. Perform massage 2-3 times a week.

Muscles of press


Fig. 18,19 - EMS electrode arrangement.

It improves blood circulation, works and tones the mouse press, promotes weight loss, lifting in problem areas of the sides and abdomen.

Facilitates the study and preparation of mice for physical activity.

Fig. 20,21- TENS electrode arrangement.

Neurostimulation covers a large area of use and penetrates deeper into body tissues.

Notice the cross-shaped arrangement of the electrodes!

Before using the device in this mode, drink 300 ml of warm water.

Install electrodes to problem places, refer to Fig.18-21. Select the massage area and type of massage "INTENSIVE," "ACUPUNCTURE" or "VIBRO MASSAGE" set the timer for 10-25 minutes. Set a comfortable intensity for you. If necessary, massage is allowed to last longer.

After using this program, it is recommended to do basic gymnastics exercises for better effect!

Muscles of legs & buttock


Figure 22 - EMS electrode arrangement - improves blood circulation, processes and tones mice, contributes to weight loss, tightening in the problem area. Facilitates the study and preparation of mice for physical activity. Used to relieve fatigue, maintain the tone of mice.

Before using the device and these programs - drink 300 ml of warm water.

Install electrodes to problem areas (see Fig.22-24), select massage area and "ACUPUNCTURE" or "VIBRO MASSAGE" modes for 10-15 minutes.

Fig. 23,24 - EMS electrode arrangement.

Facilitates the study and preparation of mice for physical activity.

Select the massage area and type of massage "ACUPUNCTURE" and set the timer for 15 minutes or "VIBRO MASSAGE" for another 10 minutes. Massage 2-3 times a week.


Figure 25.28 - EMS electrode arrangement - improves blood circulation, processes and tones joints and muscles.

Select the massage area and type of massage "ACUPUNCTURE", "INTENSE" and set the timer for 15 minutes or "VIBRO MASSAGE" for another 10 minutes. Perform massage 2-3 times a week.

Figure 26.27 - this EMS method of using electrodes, improves blood circulation, processes and tones the mice, relieves voltage.

Facilitates the study and preparation of mice for physical activity.

Select the massage area and type of massage "ACUPUNCTURE" "INTENSIVE," "VIBRO MASSAGE" or "THUMP" set the timer for 15-20 minutes.

Massage 2-3 times a week.

You can also choose from accessories - KNEE/ELBOW or SLIPPERS on any of the massage modes!

"Fitness" and "weight control" programs

This device has "FITNESS" modes to maintain the tone of mice and "WEIGHT CONTROL" for control of weight.

Before using the device and these programs - drink 300 ml of warm water.

Install electrodes to problem places (refer to Fig.18-21 and Fig.22-24) or use "WAIST" accessory.

Select "FITNESS" or "WEIGHT CONTROL" program and set the timer for 10-25 min.

Set a comfortable massage intensity for you. If necessary, massage is allowed to last longer.

Massage 2-3 times a week.


After using this program, it is recommended to do basic gymnastics exercises for better effect!

Note: to achieve the best results, these programs must be used in combination with physical activity and a balanced diet.

POWER SAVING FUNCTION

The EMS-750 device has an energy saving function.

Display is switched off during 2-3 min.

CARE AND MAINTENANCE

Device Care

Clean the device body with dry cloth only. Make sure that the device is turned off before cleaning. If salted areas appear on the LCD display, wipe it with a slightly wet soft cloth and wipe dry. Do not use chemical or abrasive cleaning and care products.

Electrode Care

Inspect the electrodes before each use. Reusable electrodes can cause mild skin irritation, lose grip, and yield less efficacy if misused. For hygienic reasons, only one person can use each pair of electrodes. To increase the life of the electrodes, before use, it is necessary to wash, degrease with an alcohol-containing liquid and dry the area of the skin on which the electrodes will be applied, as well as remove excess hair on this area of the skin. To ensure the longest possible adhesion of the electrodes, they should be carefully cleaned with a wet, non-fiber cloth or clean the lower surface of the electrodes with a jet of warm water, then shake the water very carefully and let dry in the air.

Disconnect the connection cable from the electrodes before cleaning under water. After use, stick the electrodes back to the film. At the end of the procedure, glue the protective film to the adhesive side of the electrode. The adhesive surface of the electrodes should be washed once after 30 applications or earlier, in case of their severe contamination and poor adhesion to the skin. Store the electrodes in a reusable plastic bag. Disconnect the electrodes as specified in the instruction. If one of the pair of electrodes is disconnected during the procedure, the pulse intensity will go to the lowest level, it is recommended to press the "ON/OFF" button and fix the electrodes again, then press the "ON/OFF" button and set the desired intensity. Do not pull the electrodes by the wires. This can damage both wires and electrodes.

UTILIZATION

Improper disposal of devices causes irreparable damage to the environment. Do not throw out unfair devices along with household waste. For this purpose, contact the dedicated instrument disposal point. Addresses of household appliances reception points for recycling can be obtained from municipal services of your city.

TECHNICAL CHARACTERISTICS

Model: EMS-750

Output curve shape: two-phase rectangular pulses

Power: 0.3 W

Output power: < 0.2 W

Output current: < 1 mA

Output voltage: 40V-180V

Pulse duration: from 30 to 350 μ s

Pulse frequency: from 1 to 290 Hz

Intensity adjustable, ranging from 0 to 90 mA (1000 ohms)

Power supply: 220 V, AA type batteries

Dimensions of the instrument, cm: 15 * 10 * 5 cm

Weight of myostimulator: 500 g.

Number of channels: 2

Number of intensity levels: 20

Number of electrodes: 10

Instrument material: plastic

Electrodes:

Electrode working surface dimensions, cm: 9, 5 * 6 s m

Length of connecting wires: 147 cm

Material: silicone

Slippers:

Slippers size, in cm: 26 cm

Slippers material: rubber

Waist:

Waist size: 121.5 * 15 cm

Grasp of the belt: 121 cm

Material: fabric cotton

Knee/elbow:

Size of knee/elbow: 67.5 * 19.5 cm

Knee/elbow grasp: 50 cm

Material: fabric cotton

Operating conditions: 15 ° C - 30 ° C
at relative humidity

Operating conditions: 15 ° C - 30 ° C
at relative air humidity 30 - 85%

Storage conditions: 0 ° C - 50 ° C
at relative air humidity up to 90%

Weight of full set: 1.9 kg.

FAULTS

FAULTS	REASON	SOLUTION
The device doesn't work	Incorrect polarity of power supply units installation	Install the power supplies while maintaining polarity
Turn off very quickly	Electrodes installed incorrectly	Put electrode's stickers correctly
The device turns off during the working mode	Weak or discharged power supplies	Install new batteries type AA, or connect to socket the power cord
	The massage finished because timers ends after 15 min	Try to make massage again or turn off the device
	Power cord/adaptor is broken	Change the power cord or batteries AA
The adhesive surfaces of the plates are too adhesive	Transparent plastic film not removed	Remove the film from the adhesive surface of the electronic plate
	Plate applied to skin immediately after washing	Dry the electrodes
	plate adhesion surface is damaged	Change the electrodes
Electrode stickers too much sticky	Electrodes used not on dry body	Put electrodes in the fridge for a night
	Electrodes washed too long or often	
	Electrodes were stored on high temperatures or directly under sunshine	
Feeling of stimulation effect is poor	Electrodes installed badly	Put electrodes densely
	Electrode overlaps one another	Apply electrodes to the body without overlapping
	Bad Contact, in the place where power cord connected	Check again good connection with powercord
	Low installed intensity	Increase the intensity of stimulation, using channel 1 or 2
Redness or irritation of the skin when stimulated	Dry or dirty surface of electrodes	Wash the adhesive surfaces of the plates for 3 seconds with the fingertips under a weak stream of water
	Damage of electrode's surface	Change electrode - stickers

WARRANTY CARD

<p>Warranty terms</p> <p>1. In the case of a malfunction during the Warranty Period, this product is subject to exchange by the retailer that sold this product.</p> <p>2. Warranty replacement is made in the presence of a completed warranty card and subject to the rules of operation described in the instructions for use.</p> <p>3. The product warranty does not apply in the following cases: - mechanical damage; - failure of the product due to ingestion of foreign objects and liquids, insects, etc.;; - use of the product in conditions and modes that differ from domestic ones.</p> <p>4. The warranty does not cover accessories and components.</p> <p>5. The warranty is also void if during the warranty period, the defective product was repaired by unauthorized persons.</p> <p>6. The warranty period is 1 year. 7. Shelf life - 3 years. The date of manufacture is indicated on the packaging.</p> <p>Manufacturer: Shenzhen As Tec Technology Co., Ltd Address: 8E Xin Bao Yi Industrial Bld, Hou ting Village, Bei Ting Rd, Shajing, ShenZhen, China</p> <p>Importer / Authorized Organization: LLC "SMARTTAIM" ADD: 129090 Moscow, Protopopovskiy pereulok, 25, B., floor 4, 1, room 28</p> <p>Customer service support: https://www.planta.ru.com +7 (495) 137-80-32</p> <p>The goods certified</p> 	<p>WARRANTY CARD</p> <p>Device: _____</p> <p>Model: _____</p> <p>Date of sale: _____</p> <p>Seller organisation : _____</p> <p>I received the product, I have no complaints about the configuration and appearance, I agree with the terms of the guarantee. Buyer's full name and signature: _____</p> <p>Opened the package, checked and sold. Seller's name: _____</p> <p>Seal.</p>
--	--

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